

Low Backgrounds with Noblesse

- New source materials and cleaning methods have enabled us to reduce the mass 36 background in Noblesse. This is desirable since ^{36}Ar is often used to indicate the amount of atmospheric contamination to a sample. In $^{40}\text{Ar}/^{39}\text{Ar}$ dating, this is a critical correction.
- We present data showing the H^{35}Cl interference to mass 36 as low as 2×10^{-15} ccSTP. H^{35}Cl could be measured both with the instrument pumping, and also in static mode with Ar present.
- For small samples, this low level of interference is similar to the uncertainty in the ^{36}Ar measurement; for larger samples it would be lost in the noise.

Introduction

In $^{40}\text{Ar}/^{39}\text{Ar}$ dating, the background at mass 36 can be especially important, as the small ^{36}Ar peak is used for trapped air correction. The mass 36 peak is composite, composed of ^{36}Ar , H^{35}Cl ($m/\Delta m = 3937$ relative to ^{36}Ar) and C_3 ($m/\Delta m = 1108$ relative to ^{36}Ar).

In practice C_3 is easily resolved with Noblesse, so that any potential interference to ^{36}Ar derives from H^{35}Cl only (Figure 1). Noblesse can resolve H^{35}Cl from ^{36}Ar in high-resolving power mode, albeit with some loss of sensitivity. Alternatively, the instrument can be operated at higher sensitivity, when H^{35}Cl is included in the ^{36}Ar measurement and a correction is necessary.




Figure 1: Mass 36 peak with and without Ar in the instrument (static and pumping, respectively). Scans made in high resolving power mode, with MRP ~6000 and sensitivity ~30000 atoms/cps, in order to observe the H^{35}Cl .

Lower Backgrounds

Following developments in novel source materials and the introduction of new cleaning methods, Noblesse now exhibits extremely low mass 36 backgrounds. Over a period of 6 months, the filaments were deliberately exchanged frequently, and the instrument baked out after each change. In total six filaments were tried.




Figure 2: Observed mass 36 background (H^{35}Cl) after bakeout. Different symbols indicate different filaments and bakeouts. As expected, bakeout increases the background, which then decreases with time.

The amounts of ^{36}Ar and H^{35}Cl were also measured (using the method described by Saxton, 2015), both for instrument background (pumping) and also with small amounts of Ar in the instrument.

The H^{35}Cl remains measurable in the presence of a small amount of Ar, but the amount is similar to the uncertainty in the ^{36}Ar measurement (Figure 3). For larger amounts of ^{36}Ar , the H^{35}Cl will be lost in the measurement noise (i.e. 2 cps of H^{35}Cl in comparison to the ^{36}Ar noise, which increases as $\sqrt{\text{signal}}$).

Given such low levels of H^{35}Cl , we expect that Ar measurements at lower resolving power (i.e. including H^{35}Cl but applying a small correction) should be satisfactory for almost all samples. This would also offer the advantage of using Noblesse at higher sensitivity.




Figure 3: Amount of H³⁵Cl observed with instrument pumping, and also in static mode. The H³⁵Cl remains measurable in static mode, but the amount is similar to the measurement uncertainty: The dashed line indicates the average background value of 1.6 cps, which corresponds to 2×10^{-15} ccSTP (Ar equivalent). For larger samples, H³⁵Cl would be lost in the noise.

References

J M Saxton, A method for measurement of ³⁶Ar without H³⁵Cl interference. *Chemical Geology*, 409 (2015), 112-117.