

Key elements for drug test reporting

Immunoassay

Cocaine metabolite (benzoylecgonine) screen, Urine, Qualitative

Common drug name

Matrix tested

Type of reporting

Presumptive positive

Methodology: Immunoassay

Positive cutoff: 300 ng/mL

This screening test was performed using immunoassay technology, which may produce false-positive or false-negative results. For questions, or to request confirmatory testing, contact the laboratory.

Not confirmatory

Methodology

Cutoff used

Limitations and interpretive guidance

Mass spectrometry

Opiates confirmation, Urine, Quantitative

Matrix tested

Type of reporting

Hydromorphone: 514 ng/mL

Methodology: Mass spectrometry

Positive cutoff: 20 ng/mL

Hydromorphone may be present after use of hydromorphone, or as a metabolite of hydrocodone or morphine.

[Insert lab's laboratory developed test statement here]

Methodology

Cutoff used

Limitations and interpretive guidance