

Self-Monitoring Blood Glucose Using Accu-Chek Performa Glucometer



Step 1:
Wash and dry
your hands.



Step 2:
Check the expiry date
of the test strip
indicated on the
container. Do not use
expired test strips.



Step 3:
Remove one test
strip from the
container and
close the cap
tightly.



Step 4:
Insert the test strip
into the meter until
it beeps. Wait for a
flashing blood drop
symbol to appear.




Step 5:
Press the lancing
device firmly against
the edge of your
fingertip. Push the
release button.



Step 6:
Relax your finger
after pricking.
Gently massage
until a drop of
blood appears.



Step 7:
Apply the blood
droplet to the yellow
area at the end of
the test strip. The
result will appear
after 5 seconds.



Step 8:
Press the finger with
a tissue to stop the
bleeding.





Step 9:
Slide the lever over
and back once to
replace a new
lancet for your next
blood glucose test.



Step 10:
Record your readings
in the blood glucose
logbook.



Step 11:
Press the "on" button
to display the home
screen, then  or  to scroll
through all your results.