


Self-Monitoring Blood Glucose Using Freestyle Optium Neo Glucometer




Step 1:
Wash and dry
your hands.



Step 2:
Check the expiry date
of the test strip. Do not
use expired test strips.



Step 3:
Insert the test strip and wait for
the blood drop symbol.



Step 4:
Remove the protective cap
by twisting 90 degrees from
the SF Lancet



Step 8:
The result
will appear
after 5
seconds.



Step 7:
Apply the blood drop
to the white area at the
end of the test strip.



Step 6:
Relax your finger after pricking.
Gently massage until a drop of
blood appears.



Step 5:
Place the end of the SF Lancet
on your fingertip. Gently push
against the fingertip.



Step 9:
Press the finger with a
tissue to stop the bleeding.



Step 10:
Dispose used SF Lancet
into a hard container
with a cap.



Step 11:
Record your readings in
the blood glucose logbook



Step 12:
Press to display the home screen,
then press & or to scroll
through all your results.