

CROZIER

Museum Crate



Ideal For

Museum-level protection for high-value & sensitive assets.

Design

Every build custom designed in alignment with exacting standards. Asset consolidation for sustainability.

Modification

Further customizable at the component level.

Reuse

Yes, long lifecycle.

Access

Loading oriented for best access. Options include one lid, multiple lids, or all walls removable.

Closures

Steel bolt plates for extended reuse.

Vapor Barrier Lining

Available upon request.

Interior Bracing

Lumber & MDO as needed.

Exterior Finish

Rounded edges, screw holes filled, all sanded. Paint or clear water-based lacquer exterior seal.

Exterior Marks

Lid registration, up-arrows, "fragile", piece number, unique serial number, dimensions, weight, ISPM 15 "bug stamp".

Wall Material

1/2" MDO (medium density overlay), two-sided.

Reinforcements

Fully battened with 3/4" pine lumber. Corner battens to reinforce joins.

Skids

Solid wood, 3.5" tall. Patterned for four-way handling access whenever possible.

Handles

Solid wood, consistent height for safe lifting.

Insulation

High-grade 2" thermal insulation foam all around. Lid gasket.

Cushioning

2" minimum of high-grade foam all around. Calculated to target the optimal static load case-by-case.

Interior Packaging

MDO, A-C plywood, 1/2" foam board, or 1/2" honeycomb for interior boxes and panels.

Documentation

Pack list, packing instructions as needed, carbon estimate, and fully detailed specs kept on file.

TSA Screening

Available upon request or regulation.

[Consult a Crating Expert](#)